

May 2, 2013

Protocol for Creating Multiplexed miRNA Libraries for Use in Illumina Sequencing (v3.2)

This protocol is a modification of the method described in:

“Barcoding bias in high-throughput multiplex sequencing of miRNA.” Alon S. et al. *Genome Res.* 2011 Sep;21(9):1506-11. Epub 2011 Jul 12.

Contact information

Brian Roberts
Hudson Alpha Institute for Biotechnology
601 Genome Way
Huntsville, AL 35806
Telephone: 256-327-5220
Email: broberts@hudsonalpha.org

Richard Myers
Hudson Alpha Institute for Biotechnology
601 Genome Way
Huntsville, AL 35806
Telephone: 256-327-0431
Email: rmyers@hudsonalpha.org

Overview

This protocol describes the preparation of multiplexed (barcoded) libraries of micro RNA (miRNA) from total RNA samples suitable for sequencing on the Illumina HiSeq and GAII platforms. The total RNA must be prepared by a technique that captures short RNA species (15nt-25nt). Acceptable techniques are phenol-chloroform extraction followed by ethanol precipitation or NorGen columns, amongst others.

MicroRNA species in the samples have an adaptor oligo (referred to as the 3' adaptor) ligated to their 3' ends. Next a different oligo (referred to as the 5' adaptor) is ligated to the 5' end. The 3' adaptor provides a binding site for a complementary RT primer. This allows for cDNA to be made from the miRNA-adaptor complex via reverse transcription. The cDNA is then used as a template for several rounds of PCR. The PCR primers have long tails (~30 nt) that extend the length of the product. The tails contain the barcoding sequences (with an index read primer binding site), the Illumina sequencing primer site, and the Illumina cluster-generating sequences.

Considerations

The first step, the 3' ligation, is probably the most important step in the protocol. It hinges on the use of a truncated form of T4 RNA ligase. The truncation renders the enzyme unable to use ATP for energy and instead must use an already adenylated oligo as a substrate. Consequently, the 3' adaptor oligo is adenylated on its 3' end. If a fully functional RNA

ligase capable of using ATP were used in this step, it would ligate the various RNA species present into concatamers, instead of only ligating the adaptor to the target RNA species. However, one must appreciate that all RNA species are targets, not just the miRNAs, leading to the formation of numerous unintended ligation products. Not only must these products be removed prior to sequencing, the other RNA species distract the adaptor from the miRNA population. Although it is difficult to calculate directly, the effective efficiency of ligation of the miRNAs present, in terms of the percent of miRNAs that actually get ligated to a 3' adaptor, is likely low. The 3' adaptor also has a 3-carbon spacer on its 5' end. This is to prevent RNA from being ligating to its 5' end in the subsequent 5' ligation reaction, which uses full length T4 RNA ligase.

Although the unintended ligation products that occur when an RNA molecule present in the sample other than a miRNA is ligated are somewhat problematic, the most problematic product formed in this protocol arises when un-ligated 3' adaptor is ligated to the 5' adaptor creating an adaptor dimer in the second ligation step. This creates a short product that following PCR is highly complementary to the intended miRNA ligation product. In fact, they only differ by the internal ~22 bp of the miRNA. The adaptor dimer will hybridize efficiently to the intended miRNA product, making the separation of the two difficult. This problem is dealt with in three ways. First, after the 3' ligation reaction, the entire reaction is ethanol precipitated. The precipitation removes the buffer of the 3' ligation reaction, which improves the efficiency of subsequent steps. Most importantly, the precipitation reduces the formation of the adaptor dimer. Why this occurs is unknown, but is proven experimentally. Second, after the ethanol precipitation, the RT primer is hybridized to the 3' ligation product. Since the RT primer is complementary to the entire length of the 3' adaptor, the hybridization serves to bind some of the un-ligated 3' adaptor, preventing adaptor dimer formation in the 5' ligation reaction. While both of these techniques reduce the formation of the adaptor dimer, much still persists and is present after PCR. It must be separated from the intended product by gel electrophoresis. However, because of the strong hybridization between the adaptor dimer and the intended miRNA product, the gel must be run under extremely denaturing conditions. To accomplish this, 10% acrylamide TBE-Urea gels are used. Furthermore, the gels are run in pre-heated buffer (90C). Although it is inconvenient to run the hot gels, our studies have shown that the 10% acrylamide TBE-Urea gels run at room temperature are not sufficiently denaturing for this application.

Detailed Procedure

General Notes

For those steps in which multiple components are added to the reaction, best practice is to make a "master mix" of the components sufficient for all reactions being performed. The samples throughout the course of the protocol should always be kept on ice or at ice temperature when not being otherwise incubated. In the development of this protocol, the samples were kept in a metal block that was kept cool in a refrigerator when not in use. (For simplicity, the protocol will say "on ice") Furthermore, the T4 RNA Ligase 2,

May 2, 2013

truncated; the RNase Inhibitor, murine; the T4 RNA Ligase 1, the SuperScript II, and the Phusion PCR Master Mix should all be kept on ice. In the following step, **“STOPPING POINT”** is written at points where the protocol can be stopped overnight. This protocol can be completed in 3 days. If the precipitation overnight incubations at -30C are shortened to 2-hour incubations at -80C, the protocol can be done in 2 busy days.

All of the reagents used along with vendors and catalog numbers are listed in Table 1 at the end of the document.

Sequences for the adaptors, RT primer, spike-in's, and PCR primers (suitable for sequencing on the Illumina HiSeq platform) are listed in Table 2 at the end of the document. The PCR primers used in Alon et al. (2011) are listed along with a set of 48 primers with 7-base barcodes designed by Dr. Shawn Levy and the Genomic Services Lab (GSL) at the Hudson Alpha Institute for Biotechnology.

3' Ligation

1. Prepare the following stock buffer, called “2X 3' Ligation Buffer”. This recipe is sufficient for many reactions and does not need to be prepared fresh each time the protocol is run. Store at -20C between uses:

250 uL 50% PEG 8000 (from T4 RNA Ligase 1 kit)
200 uL 10X T4 RNA Ligase Buffer (from T4 RNA Ligase 1 kit)
550 uL DNase, RNase free water

2. Make a stock solution of the spike-in controls. Make a large batch suitable for multiple runs of this protocol. Store at -80C. The concentrations listed here are suitable for the tissue types used in the development of this protocol (human brain, liver, and heart). However, it is expected that the total input of the spike-ins will need to be adjusted for different sample types.

Final Concentrations:

5000 pM miRNASeq Multiplex 22bp Spike In
500 pM miRNASeq Multiplex 25bp Spike In
50 pM miRNASeq Multiplex 30bp Spike In

3. Combine the following in a 0.2 mL PCR tube. The required input mass of RNA for the protocol is difficult to precisely state. Human brain at 200 ng has been used successfully. However, brain is a miRNA-rich tissue. It is expected that different tissues will have differing lower input limits depending on the ratio of miRNAs to total RNA in that tissue. For immortalized human cell lines (i.e. GM12878 and K562) we used 500 ng of total RNA with good results.

1 uL 10uM miRNASeq Multiplex 3' Adaptor
1 uL Spike In stock (from step 2)

May 2, 2013

200-500 ng total RNA

DNase, RNase free water to total volume of 7 uL

4. Gently mix by flicking the tube and spin down the tube in a tabletop mini-centrifuge. Incubate for 2 min at 70C in a pre-heated thermal cycler. Immediately chill on ice following incubation.
5. To each sample add the following:

10 uL 2X 3' Ligation Buffer
2 uL T4 RNA Ligase 2, truncated
1 uL RNase Inhibitor, murine
6. Gently mix the components and spin down. Incubate for 1 hour at 25C in a thermal cycler. (Note: Incubation times longer than 1 hour have been shown to produce undesired products.)
7. Ethanol precipitate the samples. Add 2.2 uL of 3M Sodium Acetate and mix well. Transfer the entire sample into a 1.7 mL microcentrifuge tube. Add 66.6 uL of 100% molecular biology grade ethanol. Mix well and spin down briefly. Incubate at -30C overnight. (Note: Longer than overnight incubation does not have adverse effects. Incubation at -80C for two hours works albeit with slightly lower yields.)
STOPPING POINT (The samples can be left in the precipitating conditions at -30C for several days.)

Resuspension, RT Primer Hybridization, and 5' Ligation

1. Remove samples from -30C incubation and gently thaw on ice. Centrifuge the samples in a refrigerated centrifuge (4C) for 20 min at ~14,000 rpm (max speed). Make sure to place the tubes in the centrifuge such that the hinges are facing away from the center of the rotor. This is important because the pellet will not be visible. By placing the hinge outward, the pellet will form near the bottom of the tube on the hinge side of the tube wall. This will allow one to more easily pipet off the supernatant without disturbing the pellet.
2. While the samples are in the centrifuge, chill an aliquot of 80% ethanol by place in it in ice water or by some other suitable method.
3. After centrifugation, pipet off the supernatant. Using a p200, place the tip of the pipet near the bottom of the tube away from the hinge side and gently remove the liquid. Add 100 uL of the chilled 80% ethanol and centrifuge again in a refrigerated centrifuge (4C) for 10 min at 14,000 rpm (max speed).

May 2, 2013

4. Again carefully remove the supernatant with the p200 as described above. After removing as much as possible with the p200, use a p20 to get the remainder, leaving behind as little liquid as possible.
5. Resuspend each sample by adding the following (a master mix can be prepared depending on the number of samples):

2 uL 10X T4 RNA Ligase Buffer (T4 RNA Ligase 1 kit)
1 uL RNase Inhibitor, murine
17 uL DNase, RNase free water
6. Transfer the resuspended samples to a clean 0.2 mL PCR tube. Add 1uL of 10 uM miRNASeq Multiplex RT Primer. Incubate as follows in a thermal cycler:

75C for 5 min
37C for 30 min
25C for 15 min
4C hold
7. While the samples are incubating, thaw the 40 uM miRNASeq Multiplex 5' Adaptor. Once thawed, incubate the adaptor at 70C for 2 min and then immediately chill on ice.
8. When the samples are finished incubating, transfer them to ice. Add the following:

0.64 uL T4 RNA Ligase 1
1 uL RNase Inhibitor, murine
0.86 uL RNase , DNase free water
1 uL 40 uM miRNASeq Multiplex 5' Adaptor
1 uL 10X T4 RNA Ligase Buffer (T4 RNA Ligase 1 kit)
1 uL 10 mM ATP (T4 RNA Ligase 1 kit)
9. Mix gently and spin down briefly. Incubate the samples for 1 hour at 25C in a thermal cycler. **STOPPING POINT** (The samples can be placed in -80C and left overnight after this step, although it is ideal to take the samples through reverse transcription before stopping)

Reverse Transcription and PCR

1. Setup the following reaction. The protocol up this point has generated ~ 26.5 uL of ligated product. Only 11 uL of the product is carried forward, so that the remainder is available for a repeat if needed. The unused product should be stored at -80C.

4 uL 5X FS Buffer (SuperScript II kit)
2 uL 0.1 M DTT (SuperScript II kit)

May 2, 2013

1 uL Deoxynucleotide Mix (10 mM each)
1 uL RNase Inhibitor, murine
1 uL SuperScript II (SuperScript II kit)
11 uL ligation product (from previous step)

2. Incubate the samples in a thermal cycler as follows:

42C for 50 min
70C for 15 min
4C hold

3. Add the following to each sample:

25 uL Phusion High-Fidelity PCR Master Mix
2.5 uL 20 uM miRNASeq Multiplex R Primer

To each individual sample add 2.5 uL of one of the twelve different indexed miRNASeq Multiplex F Primers at 20 uM, being sure to note which sample received which barcoded primer. Mix the samples and spin down.

4. Incubate the samples in a thermal cycler as follows:

94C for 30s

15 cycles of:
94C for 10s
72C for 45s

65C for 5 min
4C hold

STOPPING POINT (The samples can be stored at -20C)

Concentration, Gel Separation, and Purification

The gels run in this protocol are the Mini-PROTEAN format from BioRad and run in the Mini-PROTEAN Tetra Cell gel system. It is expected that using a different gel system would require that extensive modifications be made to this protocol.

1. Transfer each sample to a 1.7 mL microcentrifuge tube. Add 250 uL of Buffer PB (MinElute Kit). Mix well and transfer to a MinElute column placed in a 2mL collection tube. Centrifuge for 1 min at max speed. Discard flow through.
2. Add 750 uL of Buffer PE (MinElute kit, ensure ethanol has been added) to the MinElute column. Centrifuge for 1 min at max speed. Discard flow through and

May 2, 2013

place column back into the same collection tube. Centrifuge again for 1 min at max speed.

3. Transfer the column to a clean 1.7 mL microcentrifuge tube. Add 17.5 uL of RNase, DNase free water. Let stand for 5 min. Centrifuge for 1 min at max speed. Discard column, keeping the flow through in the microcentrifuge tube.
4. To each sample, add 17.5 uL of 2X TBE-Urea Sample Buffer. Mix well and spin down. Set the samples aside at room temperature.
5. Prepare DNA ladder working solutions. This recipe makes enough for several runs and need not be made fresh. Store at 4C.

20 bp Ladder

200 uL 2X TBE-Urea Sample Buffer

180 uL DNase, RNase free water

20 uL 20 bp DNA Ladder stock solution (Bayou BioLabs)

100 bp Ladder

200 uL 2X TBE-Urea Sample Buffer

190 uL DNase, RNase free water

10 uL 100 bp DNA Ladder stock solution (NEB)

6. At this point in the protocol, a hot gel will be run. Since this involves using near-boiling TBE buffer, extreme caution should be used. Additionally, protective equipment such as aprons and gloves should be worn.
7. Preheat a heating block to 95C.
8. Make 1X TBE buffer from 10X TBE buffer stock. Make 1 liter, sufficient for one or two gels. A single gel can accommodate four samples with no spacer lane between samples. Each sample will be split and run in two lanes to avoid interference from the adaptor dimer. It is not recommended to run more than two gels at a time.
9. Pre-warm 10% TBE-Urea Mini-PROTEAN gel(s) in hot tap water (no hotter than what comes out of the tap). Leave them in their packaging and weigh them down so they don't float. Also, warm the gel holder in the water.
10. In a microwave, heat 900 mL of 1X TBE buffer split into aliquots of 450 mL in two 500 mL Pyrex beakers with Saran wrap partially covering the top to 90C-95C. Heat in increments of 2-5 min (depending on microwave power). Between heating increments, carefully stir the buffer with a thermometer and check the temperature. Do not boil the buffer.

May 2, 2013

11. When the heating of the buffer is nearing completion, place the samples into the preheated heating block at 95C. Also place the 20 bp and 100 bp working solutions in the heating block. Ensure that every sample resides at 95C for at least two minutes before it is loaded onto the gel. It is not detrimental for the samples to remain in the heating block for more than 2 min, up to ~30 min.
12. Remove the gels and gel holder from the warm water. Remove the gels from their packaging, ensuring to remove the green tape at the bottom of the gel and the lane comb. Assemble the gels in the gel holder.
13. Pour the now hot 1X TBE buffer (90C-95C) into the gel assembly, filling it to the top.
14. With a p20 set to 15 uL, pipet up and down in each well of the gel. This is to remove any urea that often crystalizes in the wells during storage. Remove any bubbles in the wells.
15. Remove the two ladder tubes (carefully, they are hot). Spin them down briefly in a tabletop mini centrifuge. Add 15 uL of the 20 bp ladder to lane 1 of the gel, pipetting carefully to avoid contamination of other lanes. The tube may make a “pop” when opened. Add 15 uL of the 100 bp ladder to lane 2.
16. Remove a pair of sample tubes from the heating block. Spin them down briefly in a tabletop mini centrifuge. Load two 15 uL aliquots of each sample into two adjacent lanes of the gel. Repeat for all of the samples. Work quickly because the gel is cooling, but carefully and deliberately.
17. Once all the samples are loaded, gently place the gel assembly into the gel box. Re-heat the remaining 1X TBE buffer to 90C in the microwave. Pour all the remaining 1X TBE into the gel box (not inside the gel assembly).
18. With a 10 mL pipet, top-off the buffer inside of the gel assembly with buffer in the gel box, filling it as near to the top as possible. This is important because the hot buffer will evaporate during the course of the run.
19. Begin running the gel at 200V. Closely monitor the current. If the current begins to rise more than 10 mA from the initial current (this is likely to happen), turn the voltage down to 190V. Continue to monitor the current and adjust the voltage lower until the current stabilizes. However, do not run the gel below 160V. The current rises because the gel and buffer are hot. The conductivity of the system is much higher than when run at room temperature. The increased conductivity allows more current to flow, which in turns heats the gel, further increasing conductivity, and creating a positive feedback loop. Thus, the current must be monitored closely during the run. Under these conditions, the gel should be run for 45 minutes.


20. Turn off the power source and disassemble the gel box. Allow the gels to cool on the bench top prior to opening their plastic cases. While the gels are cooling, for each gel, add 50 mL of 1X TBE to a suitably sized gel staining container. Add 5 μ L of SYBR Gold 10,000X stock to each 50 mL TBE aliquot and mix. Wrap the container in aluminum foil to protect it from light. Open the plastic case of the now cooled gel and place the gel into the staining container with the TBE and SYBR Gold. Re-cover the container with the aluminum foil and rock on a gel rocker for 10 minutes.
21. While the gel is staining, prepare the following for each sample. With a 20-gauge needle, poke a hole in the bottom of a 0.5 mL microcentrifuge tube. Place this tube into a 1.7 mL centrifuge tube.
22. Place a sheet of Saran wrap on a UV-transilluminator. Transfer the gel from the staining solution onto the Saran wrap sheet. Capture an image of the gel under UV illumination with an appropriate gel visualization system (i.e. UVP EC3 Imaging System).
23. Transfer the gel by picking up the Saran wrap to a UV-transilluminator that can be accessed for subsequent gel excision steps (maybe the same as where the image was taken). With razor blades and forceps, carefully excise the 135 bp band for each sample (see Figure 1). Since each sample was loaded in two aliquots in adjacent lanes, cut both bands from the same sample out together. Replace the razor blades and forceps after every time they touch the gel to avoid cross-contamination. Place the gel fragments into the 0.5 mL microcentrifuge tube with the hole in the bottom.
24. Transfer the 0.5 mL microcentrifuge tubes nested in 1.7 mL microcentrifuge tubes containing the gel pieces into a microcentrifuge. Spin at max speed for 1 min. The gel fragment should be in the bottom of 1.7 mL microcentrifuge tube in small pieces. If some of the gel fragment is retained in the 0.5 mL, spin at max speed for another minute.
25. Prepare the following stock, called "Soaking Solution". This recipe makes enough for many samples as does not need to be prepared fresh every time. Store at room temperature.
 - 2 mL 5M Ammonium Acetate
 - 2 mL 1% SDS solution
 - 4 μ L 0.5M EDTA
 - 16 mL RNase, DNase free water
26. Add 300 μ L of the Soaking Solution to each sample. Incubate with agitation at 70C for 2 hours.

May 2, 2013

27. Transfer each sample (including gel pieces) to a Spin-X Centrifuge Tube Filter, 0.22 μ m Cellulose Acetate, sitting in its accompanying microcentrifuge tube. Spin in a microcentrifuge at max speed for 1 min.
28. Transfer the flow-through to a new 1.7 mL microcentrifuge tube. Add 1 μ L of 10 μ g/ μ L glycogen. Add 300 μ L of 100% isopropanol. Vortex and spin down briefly. Incubate overnight at -30C. **STOPPING POINT** (The samples can be kept in the precipitating conditions at -30C for several days.)
29. Spin the samples in a refrigerated centrifuge (4C) for 20 min at 14,000 rpm (max speed). Again, place the hinges of the tubes outward so that the location of the pellet is predictable.
30. While the samples are in the centrifuge, chill an aliquot of 80% ethanol by place in it in ice water or by some other suitable method.
31. After centrifugation, pipet off the supernatant. Using a p200, place the tip of the pipet near the bottom of the tube away from the hinge side and gently remove the liquid. Add 100 μ L of the chilled 80% ethanol and centrifuge again in a refrigerated centrifuge (4C) for 10 min at 14,000 rpm (max speed).
32. Again carefully remove the supernatant with the p200 as described above. After removing as much as possible with the p200, use a p20 to get the remainder, leaving behind as little liquid as possible.
33. Resuspend the pellet in 10 μ L of EB buffer (MinElute Kit). Measure the concentration of the sample with a suitable method (QBit HS DNA is preferred with 1 μ L of sample input). The sample is ready for sequencing. Typically, this protocol yields 10 μ L of 1-4 ng/ μ L product, depending on sample input mass and sample type. Although it depends on the level of multiplexing, 0.5 ng/ μ L or higher libraries are concentrated enough for sequencing. If a lower yield is expected, the pellet can be resuspended in a lower volume to yield a higher concentration product.

May 2, 2013

Figure 1 – Sample Gel Image



May 2, 2013

Table 1 - Materials

Product Name	Component	Manufacturer	Cat #
T4 RNA Ligase 2, truncated	T4 RNA Ligase 2, truncated	NEB	M0242L
RNAse Inhibitor, Murine	RNAse Inhibitor, Murine	NEB	M0314S
T4 RNA Ligase 1	T4 RNA Ligase 1	NEB	M0204S
T4 RNA Ligase 1	10X T4 RNA Ligase Reaction Buffer	NEB	M0204S
T4 RNA Ligase 1	10 mM ATP	NEB	M0204S
T4 RNA Ligase 1	50% PEG 8000	NEB	M0204S
SuperScript II Reverse Transcriptase	SuperScript II RT	Invitrogen	18064-071
SuperScript II Reverse Transcriptase	5X FS Buffer	Invitrogen	18064-071
SuperScript II Reverse Transcriptase	0.1 M DTT	Invitrogen	18064-071
Deoxynucleotide Solution Mix	Deoxynucleotide Solution Mix	NEB	N0447S
Phusion High-Fidelity PCR Master Mix with HF Buffer	Phusion High-Fidelity PCR Master Mix with HF Buffer	NEB	M0531S
Ethanol	Ethanol	Sigma-Aldrich	E7023-500ML
MinElute PCR Purification Kit	MinElute Spin Column	Qiagen	28004
MinElute PCR Purification Kit	Buffer PB	Qiagen	28004
MinElute PCR Purification Kit	Buffer PE (concentrate)	Qiagen	28004
MinElute PCR Purification Kit	Buffer EB	Qiagen	28004
Sodium Acetate, 3M, molecular biology grade	Sodium Acetate, 3M, molecular biology grade	Usb	75897
100 bp DNA Ladder	100 bp DNA Ladder	NEB	N3231L
20 bp DNA Ladder	20 bp DNA Ladder	Bayou BioLabs	L-100
TBE-Urea Sample Buffer	TBE-Urea Sample Buffer	BioRad	161-0768
10% Mini-PROTEAN TBE-Urea Precast Gel	10% Mini-PROTEAN TBE-Urea Precast Gel	BioRad	456-6033
Accugene 10X TBE Buffer	Accugene 10X TBE Buffer	Lonza	50843
SYBR Gold 10,000X	SYBR Gold 10,000X	Invitrogen	S11494
5M Ammonium Acetate	5M Ammonium Acetate	Ambion	AM9070G
0.5M EDTA	0.5M EDTA	Ambion	AM9261
Sodium Dodecyl Sulfate (SDS)	Sodium Dodecyl Sulfate (SDS)	Sigma-Aldrich	L6026-50G
Spin-X Centrifuge Tube Filter	Spin-X Centrifuge Tube Filter	Costar	8161
Glycogen 20ug/uL	Glycogen 20ug/uL	Invitrogen	10814-010
Isopropanol, Molecular Biology Grade	Isopropanol, Molecular Biology Grade	Fisher	BP2618-500

May 2, 2013

Table 2 - Oligo List

Name	Sequence
miRNASeq Multiplex 3' Adaptor	/5rApp/ACGGGCTAATATTTATCGGTGG/3SpC3/
miRNASeq Multiplex 5' Adaptor	rUrCrCrCrUrArCrArCrGrArCrGrCrUrCrUrUrCrCrGrArUrCrUrC
miRNASeq Multiplex RT primer	GCTCCACCGATAAATATTAGCCCGT
miRNASeq Multiplex 22bp Spike In	rArGrCrGrCrUrUrGrCrArGrArGrArGrArArUrCrArG
miRNASeq Multiplex 25bp Spike In	rGrCrGrUrGrGrArCrArCrArUrCrUrGrUrCrGrGrCrCrArUrArC
miRNASeq Multiplex 30bp Spike In	rArArCrCrGrCrArCrArCrCrUrGrCrCrGrArUrGrUrCrUrUrCrArCrCrG

PCR primers from Alon et al. (2011)

PCR Primers	Sequence	Barcode on Primer	Barcode as Sequenced
miRNASeq Multiplex R Primer	AATGATACGGCGACCAACCGAGATCTACACTCTTCCCTACACGACGCTCTTCCGATCT		
miRNASeq Multiplex F Primer Index 1	CAAGCAGAAGACGGCATAACGAGATCGTGATGCTCCACCGATAAATATTAGCCCGT	CGTGAT	ATCACG
miRNASeq Multiplex F Primer Index 2	CAAGCAGAAGACGGCATAACGAGATACATCGGCTCCACCGATAAATATTAGCCCGT	ACATCG	CGATGT
miRNASeq Multiplex F Primer Index 3	CAAGCAGAAGACGGCATAACGAGATGCCTAAGCTCCACCGATAAATATTAGCCCGT	GCCTAA	TTAGGC
miRNASeq Multiplex F Primer Index 4	CAAGCAGAAGACGGCATAACGAGATTGGTCAGCTCCACCGATAAATATTAGCCCGT	TGGTCA	TGACCA
miRNASeq Multiplex F Primer Index 5	CAAGCAGAAGACGGCATAACGAGATCACTGTGCTCCACCGATAAATATTAGCCCGT	CACTGT	ACAGTG
miRNASeq Multiplex F Primer Index 6	CAAGCAGAAGACGGCATAACGAGATATTGGCGCTCCACCGATAAATATTAGCCCGT	ATTGGC	GCCAAT
miRNASeq Multiplex F Primer Index 7	CAAGCAGAAGACGGCATAACGAGATGATCTGGCTCCACCGATAAATATTAGCCCGT	GATCTG	CAGATC
miRNASeq Multiplex F Primer Index 8	CAAGCAGAAGACGGCATAACGAGATTCAAGTGCTCCACCGATAAATATTAGCCCGT	TCAAGT	ACTTGA
miRNASeq Multiplex F Primer Index 9	CAAGCAGAAGACGGCATAACGAGATCTGATCGCTCCACCGATAAATATTAGCCCGT	CTGATC	GATCAG
miRNASeq Multiplex F Primer Index 10	CAAGCAGAAGACGGCATAACGAGATAAGCTAGCTCCACCGATAAATATTAGCCCGT	AAGCTA	TAGCTT
miRNASeq Multiplex F Primer Index 11	CAAGCAGAAGACGGCATAACGAGATGTAGCCGCTCCACCGATAAATATTAGCCCGT	GTAGCC	GGCTAC
miRNASeq Multiplex F Primer Index 12	CAAGCAGAAGACGGCATAACGAGATTACAAGGCTCCACCGATAAATATTAGCCCGT	TACAAG	CTTGTA

May 2, 2013

microRNASeq primers with barcodes designed by Dr. Shawn Levy and the Genomic Services Lab (GSL) at the Hudson Alpha Institute for Biotechnology

PCR Primers	Sequence	Barcode on Primer	Barcode as Sequenced
GSL_miRNASeq F Primer Index 01	CAAGCAGAAGACGGCATACGAGATCAAACGCGCTCCACCGATAAAATATTAGCCCGT	CAAACGC	GCGTTTG
GSL_miRNASeq F Primer Index 02	CAAGCAGAAGACGGCATACGAGATATGAGCCGCTCCACCGATAAAATATTAGCCCGT	ATGAGCC	GGCTCAT
GSL_miRNASeq F Primer Index 03	CAAGCAGAAGACGGCATACGAGATGAATCCGGCTCCACCGATAAAATATTAGCCCGT	GAATCCG	CGGATTC
GSL_miRNASeq F Primer Index 04	CAAGCAGAAGACGGCATACGAGATGGCGTTAGCTCCACCGATAAAATATTAGCCCGT	GGCGTTA	TAACGCC
GSL_miRNASeq F Primer Index 05	CAAGCAGAAGACGGCATACGAGATGGCTCAAGCTCCACCGATAAAATATTAGCCCGT	GGCTCAA	TTGAGCC
GSL_miRNASeq F Primer Index 06	CAAGCAGAAGACGGCATACGAGATTACGCTGCTCCACCGATAAAATATTAGCCCGT	TCACGCT	AGCGTGA
GSL_miRNASeq F Primer Index 07	CAAGCAGAAGACGGCATACGAGATCTGCGTTGCTCCACCGATAAAATATTAGCCCGT	CTGCGTT	AACGCAG
GSL_miRNASeq F Primer Index 08	CAAGCAGAAGACGGCATACGAGATGCTTAGGGCTCCACCGATAAAATATTAGCCCGT	GCTTAGG	CCTAAGC
GSL_miRNASeq F Primer Index 09	CAAGCAGAAGACGGCATACGAGATAGGGTAGGCTCCACCGATAAAATATTAGCCCGT	AGGGTAG	CTACCCT
GSL_miRNASeq F Primer Index 10	CAAGCAGAAGACGGCATACGAGATTACGAGCGCTCCACCGATAAAATATTAGCCCGT	TACGAGC	GCTCGTA
GSL_miRNASeq F Primer Index 11	CAAGCAGAAGACGGCATACGAGATTCTGCGTGCTCCACCGATAAAATATTAGCCCGT	TCTGCGT	ACGCAGA
GSL_miRNASeq F Primer Index 12	CAAGCAGAAGACGGCATACGAGATCTTACCGCTCCACCGATAAAATATTAGCCCGT	CTTACCG	CGGTAAG
GSL_miRNASeq F Primer Index 13	CAAGCAGAAGACGGCATACGAGATGGGAGTAGCTCCACCGATAAAATATTAGCCCGT	GGGAGTA	TACTCCC
GSL_miRNASeq F Primer Index 14	CAAGCAGAAGACGGCATACGAGATACTTGGGGCTCCACCGATAAAATATTAGCCCGT	ACTTGGG	CCCAAGT
GSL_miRNASeq F Primer Index 15	CAAGCAGAAGACGGCATACGAGATGGGTTTGGCTCCACCGATAAAATATTAGCCCGT	GGGTTTG	CAAACCC
GSL_miRNASeq F Primer Index 16	CAAGCAGAAGACGGCATACGAGATTGAGTCCGCTCCACCGATAAAATATTAGCCCGT	TCAGTCC	GGACTGA
GSL_miRNASeq F Primer Index 17	CAAGCAGAAGACGGCATACGAGATTCGAGGTGCTCCACCGATAAAATATTAGCCCGT	TCGAGGT	ACCTCGA
GSL_miRNASeq F Primer Index 18	CAAGCAGAAGACGGCATACGAGATAATCGCGCTCCACCGATAAAATATTAGCCCGT	AATCGCG	CGCGATT
GSL_miRNASeq F Primer Index 19	CAAGCAGAAGACGGCATACGAGATCTGCGATGCTCCACCGATAAAATATTAGCCCGT	CTCGCAT	ATGCGAG
GSL_miRNASeq F Primer Index 20	CAAGCAGAAGACGGCATACGAGATCACGTAGGCTCCACCGATAAAATATTAGCCCGT	CACGTAG	CTACGTG
GSL_miRNASeq F Primer Index 21	CAAGCAGAAGACGGCATACGAGATAGCGACTGCTCCACCGATAAAATATTAGCCCGT	AGCGACT	AGTCGCT
GSL_miRNASeq F Primer Index 22	CAAGCAGAAGACGGCATACGAGATTGACCTCGCTCCACCGATAAAATATTAGCCCGT	TGACCTC	GAGGTCA
GSL_miRNASeq F Primer Index 23	CAAGCAGAAGACGGCATACGAGATCACCGAAGCTCCACCGATAAAATATTAGCCCGT	CACCGAA	TTCGGTG
GSL_miRNASeq F Primer Index 24	CAAGCAGAAGACGGCATACGAGATATGCAGCGCTCCACCGATAAAATATTAGCCCGT	ATGCAGC	GCTGCAT
GSL_miRNASeq F Primer Index 25	CAAGCAGAAGACGGCATACGAGATCAGTGGTGCTCCACCGATAAAATATTAGCCCGT	CAGTGGT	ACCACTG
GSL_miRNASeq F Primer Index 26	CAAGCAGAAGACGGCATACGAGATGCTTTCGCTCCACCGATAAAATATTAGCCCGT	GTCTTGC	GCAAGAC
GSL_miRNASeq F Primer Index 27	CAAGCAGAAGACGGCATACGAGATTGGAAGGGCTCCACCGATAAAATATTAGCCCGT	TGGAAGG	CCTTCCA
GSL_miRNASeq F Primer Index 28	CAAGCAGAAGACGGCATACGAGATGGGATCTGCTCCACCGATAAAATATTAGCCCGT	GGGATCT	AGATCCC
GSL_miRNASeq F Primer Index 29	CAAGCAGAAGACGGCATACGAGATATACGGCGCTCCACCGATAAAATATTAGCCCGT	ATACGGC	GCCGTAT
GSL_miRNASeq F Primer Index 30	CAAGCAGAAGACGGCATACGAGATTCGGATGGCTCCACCGATAAAATATTAGCCCGT	TCGGATG	CATCCGA
GSL_miRNASeq F Primer Index 31	CAAGCAGAAGACGGCATACGAGATGCTTCACGCTCCACCGATAAAATATTAGCCCGT	GCTTCAC	GTGAAGC
GSL_miRNASeq F Primer Index 32	CAAGCAGAAGACGGCATACGAGATCCTACTGGCTCCACCGATAAAATATTAGCCCGT	CCTACTG	CAGTAGG
GSL_miRNASeq F Primer Index 33	CAAGCAGAAGACGGCATACGAGATACCACGAGCTCCACCGATAAAATATTAGCCCGT	ACCACGA	TCGTGGT
GSL_miRNASeq F Primer Index 34	CAAGCAGAAGACGGCATACGAGATGGTGCTAGCTCCACCGATAAAATATTAGCCCGT	GGTGCTA	TAGCACC
GSL_miRNASeq F Primer Index 35	CAAGCAGAAGACGGCATACGAGATCAAGGCTGCTCCACCGATAAAATATTAGCCCGT	CAAGGCT	AGCCTTG
GSL_miRNASeq F Primer Index 36	CAAGCAGAAGACGGCATACGAGATAGGCGAAGCTCCACCGATAAAATATTAGCCCGT	AGGCGAA	TTCGCTT
GSL_miRNASeq F Primer Index 37	CAAGCAGAAGACGGCATACGAGATTCAGGTGGCTCCACCGATAAAATATTAGCCCGT	TCAGGTG	CACCTGA

May 2, 2013

GSL_miRNASeq F Primer Index 38	CAAGCAGAAGACGGCATAACGAGATCGTCCATGCTCCACCGATAAATATTAGCCCGT	CGTCCAT	ATGGACG
GSL_miRNASeq F Primer Index 39	CAAGCAGAAGACGGCATAACGAGATCTCGAGAGCTCCACCGATAAATATTAGCCCGT	CTCGAGA	TCTCGAG
GSL_miRNASeq F Primer Index 40	CAAGCAGAAGACGGCATAACGAGATAGTCGTGGCTCCACCGATAAATATTAGCCCGT	AGTCGTG	CACGACT
GSL_miRNASeq F Primer Index 41	CAAGCAGAAGACGGCATAACGAGATGAGCAGAGCTCCACCGATAAATATTAGCCCGT	GAGCAGA	TCTGCTC
GSL_miRNASeq F Primer Index 42	CAAGCAGAAGACGGCATAACGAGATAGATGCCGCTCCACCGATAAATATTAGCCCGT	AGATGCC	GGCATCT
GSL_miRNASeq F Primer Index 43	CAAGCAGAAGACGGCATAACGAGATACCTCACGCTCCACCGATAAATATTAGCCCGT	ACCTCAC	GTGAGGT
GSL_miRNASeq F Primer Index 44	CAAGCAGAAGACGGCATAACGAGATGATGCCAGCTCCACCGATAAATATTAGCCCGT	GATGCCA	TGGCATC
GSL_miRNASeq F Primer Index 45	CAAGCAGAAGACGGCATAACGAGATCCCTTCTGCTCCACCGATAAATATTAGCCCGT	CCCTTCT	AGAAGGG
GSL_miRNASeq F Primer Index 46	CAAGCAGAAGACGGCATAACGAGATTCTGGGCTCCACCGATAAATATTAGCCCGT	TTCCTGG	CCAGGAA
GSL_miRNASeq F Primer Index 47	CAAGCAGAAGACGGCATAACGAGATTGGACACGCTCCACCGATAAATATTAGCCCGT	TGGACAC	GTGTCCA
GSL_miRNASeq F Primer Index 48	CAAGCAGAAGACGGCATAACGAGATTGAAGGGGCTCCACCGATAAATATTAGCCCGT	TGAAGGG	CCCTTCA