

Certificate of Analysis



Human Bladder Microvascular Endothelial Cells (HBdMEC)

Description

Product Name	HBdMEC-c / HBdMEC-p
Order Number	C-12291 / C-12292
Lot Number	7032704.11
Donor Age / Sex / Race	46 / male / caucasian
Tissue / Localization	bladder / tunica muscularis
Number of Viable Cells	500.000
Freezing Medium	Cryo-SFM (Order No.: C-29910)
QC Evaluation Medium	Endothelial Cell Growth Medium MV (Order No.: C-22020)
Stage of Culture	HBdMEC-c: thawing and seeding results in passage 2 (3 rd culture) HBdMEC-p: shipped in passage 2 (3 rd culture)

Viability & Growth Characteristics

Viability	97 % * (adherence of cells could vary under customer culture conditions)
Population Doubling (PD) Time in Log Phase	35.6 h / PD *
Population Doublings	> 15 *

Phenotypic Characterization (tested within the first two passages)

CD31	positive by immunohistochemical staining
vWF / Factor VIII-related antigen	positive by immunohistochemical staining
Dil-Ac-LDL uptake	positive by immunohistochemical staining
Smooth muscle α -actin	negative by immunohistochemical staining

Test for Microbial Contaminants and Infectious Viruses

Bacteria / Fungi / Mycoplasma	negative / negative / negative
HIV-1	negative
HBV / HCV	negative / negative

* Using PromoCell's standardized culture system and procedures. The stated values may vary under customer culture conditions.

The tissue used by PromoCell for the isolation of human cell cultures is derived from donors who have signed an informed consent form, which outlines in detail the purpose of the donation and the procedure for processing the tissue (www.promocell.com/ethics).



Silke Busch, Ph.D.
Quality Control Manager

Date: 24 Sep 2008

FOR IN VITRO USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC PROCEDURES.

Orders and Technical Information

PromoCell GmbH
Sickingenstr. 63 / 65
69126 Heidelberg
Germany

Germany
France
United Kingdom
USA
International

Phone: 0800 – 776 66 23
Phone: 0800 90 93 32
Phone: 0800 96 03 33
Phone: 1 – 866 – 251 – 2860
Phone: +49 6221 – 649 34 0

Email: info@promocell.com
www.promocell.com